

# EasyDent June Newsflash



August 2016

## Enhancements


### New Perio Charting

There a brand new type of General Perio Charting available! You can watch a 10 minute Video about this at: [www.ezdent.com/Chartingvid.htm](http://www.ezdent.com/Chartingvid.htm) select the topic [General Perio Charting - Newest Version](#).

When you start using Perio Charting you will find 2 options: **Perio General Newest**, and Perio General Legacy. To use the new program, simply select the **Perio General Newest** option.

It has many new features:

- Separate measurements for each pocket, instead of for each side of the tooth.
- You may now specify Recession, Mobility, Furcation, & Bleeding
- Missing teeth and/or implants are represented graphically
- View Graph on-screen of recessions and pocket depths
- Adding text comments
- Automatically converts legacy perio charts, if applicable



# EasyDent June Newsflash

## EDR - Electronic Dental Records

- You can now create a **Second EDR Visit Note** for the **same date**. The Notes Date will automatically have a **"B"** as its suffix
- You can limit Users to only be able to use the Tablet Style Menu. From Advanced Utilities, Define User Security, check the option called **"Limit User to Tablet Style Menu"**. You might wish to set this for chairside providers that don't need the desktop menu.

## Prescriptions

- On the Prescriptions screen there is a new button **"Repeat Last Script"** you can use to quickly create a duplicate of the last script you did for a patient.
- When writing Prescriptions you can have 2 lines of Instructions per Medication, instead of just 1 line.

## Scheduling

- There is a new option you can set so that when you are on the Schedule in **"View Only"** Mode the patient's name **will not** be displayed. To set this from the Schedule, select the top Toolbar **Utilities** then **Utilities Option Menu**, then **Set Scheduling Option**, click the Option labeled **"Don't Show Patient Name in View Only Mode"**.


## Reports

- The Monthly Report **"Yearly Analysis"** has a new **YearToDate** column.
- There is new report to let you easily compare overall practice statistics of where you are year-to-date with the same period last year. From the Reports Menu select Management Reports, then **"Year to Date Practice Pulse"**.
- On the Daily Reports Menu, there is a new type of Bank Deposit called **"Bank Deposit Single Column Format"**. Please try it out.

## Enhancements List

Find out all the latest enhancements that are available by clicking your Primary Menu **"Enhance"** Light Bulb Icon and selecting **"View Available Internet Updates"** button.

[www.ezdent.com/Enhancements.htm](http://www.ezdent.com/Enhancements.htm)



Thank you for continuing to send us great ideas and your continued support.

Sincerely,

**Your EasyDent Staff** 😊