



Updating your Tablet/PC Programs from the Internet

Typically when you run the “**Download / Upload**” option and select “**Download**” you will pick up the latest program updates from the Cloud. If you ever need to get the latest program updates directly from the EasyDent web site to your local computer; then follow the steps below.

- 1) From the Cloud Menu select the top toolbar “**Tools**” option, then “**Download Internet Updates**”.




- 2) Launch the file you downloaded, “EZDENTUP.EXE” and click the “Yes” prompt.




Updating your Tablet/PC Programs from the Internet

3) Change the “X:\” folder to “C:\” then click the “Unzip” button.




4) Type in the password “YYYYup” and click the “OK” button. The YYYY is the current year like “2016up” or “2017up”, etc. This password changes generally it is the year followed by Lowercase “up”. If you have trouble with the password call Data Tec at (636) 256-7401 during business hours and ask for the update password.



Updating your Tablet/PC Programs from the Internet

- 5) When the unzip starts running you will see a progress bar. The update usually completes within just a couple minutes.



- 6) When the unzip is complete, click the "OK" button, and then close the Unzip Window.

