


# Charting

## 28.1 Charting- Restorative

From the patient's General screen, select the top toolbar "Window" option, then Charting, then "Restorative".




To make changes to the Chart click the update button. You can right click on the teeth diagrams to insert predefined images, or use the drawing tools on the left, by clicking on a drawing tool and then using your mouse or the pen of a pen enabled device, like a Tablet PC.

# Charting

## 28.2 Charting- Perio General

From the patient's General screen, select the top toolbar "Window" option, then Charting, then "Perio-General". On this chart you can click the top toolbar option "Quick Entry" to display a window that lets you enter values quickly.



Periodontal Chart Acct:00002 LINCOLN COKEY DOB: 11/18/1982

New Clear Quick Entry Print Save | Cancel Undo | Exit

C																	
B																	
A	233	456	777	8910	1011	544	444	444	455	666	677	777	212	102	212	332	1
	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	2
																	3
																	4
	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	5
A	341	221	120	122	011	221	221	121	131	132	232	112	553	141	531	422	6
B																	7
C																	8

Date: 07162010 Periodontal Record: First Complete Charting for Patient Acct:00002 LINCOLN COKEY DOB: 11/18/1982 Bleeding Sites Teeth

C																	
B																	
A	224	213	314	213	212	212	324	213	321	021	223	133	245	222	232	301	32
	32	31	30	29	28	27	26	25	24	23	22	21	20	19	18	17	31
																	30
																	29
	32	31	30	29	28	27	26	25	24	23	22	21	20	19	18	17	28
A	112	112	111	012	321	121	113	333	121	422	323	232	434	031	300	432	27
B																	26
																	25
																	24
																	23
																	22
																	21
																	20
																	19
																	18
																	17

# Charting

## 28.3 Charting- Perio Detailed

From the patient's General screen, select the top toolbar "Window" option, then Charting, then "Perio-Detailed". Click on the teeth you wish to enter info for and then click on "Facial" to change to "Lingual" and vice a versa.

EasyDent Perio Charting

2      LINCOLN      COKEY

Acct No      Patient: Last Name      First Name      Initial

Right      Left

1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16

32 31 30 29 28 27 26 25 24 23 22 21 20 19 18 17

**Facial**  
Click to Change View

	D		M	
Pocket	2	3	1	
Bleeding		B		
Recession	1	2	3	
Furcation		F		
Mobility		½		

Suppuration	S		
FGM	+1	+2	+3
MGJ	3	2	3
Attachment		2	

Missing

View Quad      Next Tooth F8

Save      Prev Tooth F7

Print

**HELP**

Cancel

**Tooth 1**

Notice you can click "View Quad" to see a quadrant at a time.