

# *EasyDent* CD Update Instructions

- If you have a newer computer with Windows XP or 2000 you may be prompted to Overwrite Files. Reply "A" for all each time you are prompted.

## Update Checklist

- 1 \_\_ Back up your computer ----- **THIS IS *EXTREMELY IMPORTANT***.
- 2 \_\_ Exit **all** running programs and/or restart your computer.  
  
**EasyDent, EasySchedule, and/or NetCom CAN NOT** be active, they must be closed!
- 3 \_\_ Insert the Update CD (**label up**) into your computer.
- 4 \_\_ Within approximately 1 minute you should be prompted with a **Menu**.  
If not then do the following:
  - Open you're My Computer icon.
  - Locate the EasyDent update CD; double-click on it, to open it.
  - Next, double-click on the **Runnow** file.
- 5 \_\_ Select the update option for **your computer**.
  - If you are at the Main Computer use the Main Computer option.  
(For dedicated servers, you can run Main Computer option at a Work Station)
  - If you are at a Work Station use the Work Station option.
- 6 \_\_ If you see an "**Abort, Retry, Fail**" message, then restart your computer, don't start EasyDent, and use the CD again, please. Another option is to reply "F".
- 7 \_\_ Remove the Update CD when the update is complete.
- 8 \_\_ Restart **EasyDent** select the Primary Menu "Enhancements Lists".

This must be accomplished at ALL Computers

### If you don't have a CD drive available on a work station you can:

- Login to EasyDent at the Work Station, from the Primary Menu; use either the top toolbar "Tools" option, then "Work Station Update" the password is SECRET, or move the mouse to the lower left hand side of Menu and select the Blue "Statistics", then pick the Red Work Station Update button, and follow all the on screen prompts.