

Warning: If you are using a Cloud or Online Backup ...

More and more offices are starting to use Cloud or Online Backup services. While these are good supplements to your local backup procedures, there are two major dangers to avoid.

1) DO NOT let Cloud or Online backup services run during the day. Running the service during the day **WILL KILL** your system performance. EasyDent has lots of large databases that are constantly changing while you use the system. Many cloud backup services default to backing up any files that change, when they change. Most systems can be set to backup specific folders or files and backup them up at scheduled times. You must schedule these types of backups to run at night.

2) Cloud or Online services are only an extra safeguard they DO NOT replace your own local backups. If your entire system crashes, (computer dies, disk crashed, etc.) the Online or Cloud Backup service can take a long time to restore from plus you don't want to be totally dependent on their service!

We recommend using www.mozypro.com and setting it to backup all folders that start with **EZW** on your main computers shared drive. Also making sure you set the backup to run overnight.

For additional performance considerations visit our website www.ezdent.com select the "Support" tab, next select "Technical notes", then select Note #32.

A direct link is:

<http://www.ezdent.com/technote/EasyDentPerformance.pdf>