




Activate NewCrop Chromium Interface

- 1) Go to **ezdent** website Support Tab, then Downloads button
- 2) Download and Install link named: [Browser CefSharp Install](#) to Main Computer (for example: X:\ changed to H:\)
- 3) Download and install Latest Updates to Main Computer
- 4) At each workstation run the workstation update **twice**
- 5) Login to EasyDent, from the Primary Menu use the “Options” button



- 6) Check the “Use Browser CefSharp eScripts Interface”.
When prompted enter the Special Password, then press “Save/Exit”

