




Activate NewCrop **WebView2** Interface

Following steps only once from any workstation

- 1) Go to EasyDent website Support Tab, then Downloads button
- 2) Download and install Latest Updates to Main Computer

Following steps must be done at each workstation that will use eScripts


- 1) Go to EasyDent website Support Tab, then Downloads button
- 2) Select Link: [Microsoft WebView2 Install](#)
- 3) Click and run the “Download” option under “Evergreen Bootstrapper”




- 4) Select Link: [Webview2 Run Times Install](#)
Download and run to Drive C:
- 5) Run the EasyDent workstation update **twice**

Following steps only once from any workstation


1) Login to EasyDent, from the Primary Menu use the “Options” button



2) Check the “Use Browser CefSharp eScripts Interface”. Note: Even though you are setting up WebView2 Browser activate the CefSharp option! When prompted enter the Standard Special Password, then press “Save/Exit”




3) From any patient’s EDR Visit Note screen use the top toolbar “eScripts” option.



continued ...

4) On the “EdgeScripts” screen select “Tools”




The screenshot shows the eScripts application window with the title bar 'eScripts: Acct# 00013'. The 'Tools' menu is highlighted with a red arrow. Below the menu, there are fields for 'User ID' (containing 'HALG'), 'Password' (containing 'XXXX'), and 'Name' (containing 'AUSTEN JANE W'). A table displays patient information:

Age	Weight Lbs	Weight Date	Height Inches	Height Date	
72	6	120	12312020	65	12312020

Below the table is a grid of colored buttons: 'New Script' (red), 'Insurance Detail & Formulary' (yellow), 'Status or Admin Pages' (pink), 'Meds Prescribed by other Doctors' (cyan), 'Renewals' (green), 'Educational Materials' (magenta), and 'Exit' (blue).

5) On the drop down, list check the option “WebView2 Mode”, the password will be the the letter “E” plus the “Special Password”.



The screenshot shows the eScripts application window with the 'Tools' menu open. The 'WebView2 Mode' option is checked and highlighted with a red arrow. The 'EasyDent' button is visible in the top right. The 'User ID' field contains 'HALG' and the 'Password' field contains 'XXXX'. The 'Name' field contains 'AUSTEN JANE W'. The patient information table is the same as in the previous screenshot:

Age	Weight Lbs	Weight Date	Height Inches	Height Date	
72	6	120	12312020	65	12312020

The grid of colored buttons is also present: 'New Script' (red), 'Insurance Detail & Formulary' (yellow), 'Status or Admin Pages' (pink), 'Meds Prescribed by other Doctors' (cyan), 'Renewals' (green), 'Educational Materials' (magenta), and 'Exit' (blue).