



Data Tec, Inc., P.O. Box 31576, Des Peres, MO 63131  
Voice (636) 256-7401 Fax (636) 230-6801

---

## Digital X-Ray & Oral Cameras

1) **EasyDent** has a direct link from patient screens to the most popular X-Ray and Imaging systems, including: **Apteryx, Dexis, EVASoft, Gendex, Kodak, Patterson, Schick, Suni, Sopro, Tiger View, XDR Radiology, Visix, and VixWin.**

2) Typically, we can program a **Free Bridge** between EasyDent and most digital x-ray & oral Camera software systems. Generally, vendors will provide specifications to program the bridge with their imaging system.

3) **DataGrabber** is a Windows' System Tray application that understands how to bridge patients that you have open in one application to another application, starting it if necessary. This provides a bridge between **EasyDent** and many popular digital x-ray systems. Go to <https://ledapteryx.com> then select products, imaging, and then DataGrabber.

4) You can quickly and easily import any windows digital image file directly into the patients **EasyDent** Electronic Dental Record. As long as the image is a standard Windows format, or the imaging vendor software provides a device driver to view their file type.

### 5) Imaging Packages

No matter what type of x-ray or oral camera system you are using, as mentioned, above we provide a direct link to two of the most popular digital imaging packages, Apteryx and Tiger View. Find them at: <https://ledapteryx.com> then select products, imaging for Apteryx and <https://www.televeresystems.com/> for Tiger View.

### 6) Scan-X

You can use Apteryx Imaging, or Patterson Imaging, Apteryx Scan-X Twain Scanner Driver, or DataGrabber. Please talk to us before selecting which approach to use.