

EasyDent

Updates via the Internet

WARNING – We always highly recommend making sure you back up before applying any changes.

1. Please Contact Data Tec and obtain the most recent unzip Internet update Password.

From a Work Station or the Main Computer.

2. Open EasyDent, from the Primary Menu select the Top toolbar “Tools” option then “Internet - Download Latest Updates”.

You will be prompted with a message telling which drive letter you will Need to use when you unzip the download file, make note of it.

3. This should link to our update file and allow you to download it.

4. Once downloaded, change the folder name from **X:** to **Z:** where **Z** is the network drive letter you use to network. (Typically H:\ or G:\, whatever letter you made note of in step 1 above.)

5. Next click the Unzip button.

6. When prompted for a password enter the unzip password you obtained from Data Tec in the first step.

7. After the files are unzipped your main computer will be updated, you **must update each work station**. Use the top Toolbar Menu “Tools” option, then “Work Station Update” to update each work station.