





eScripts using Allscripts

You can use Allscripts with EasyDent **without** having to **re-type** in any patient demographic information. www.allscripts.com/products-services/products/eprescribe



1) From the Patient's General Screen select the “**eScript**” Icon, or the “**Script**” Icon on a patient's EDR Visit Note screen.



2) The Allscripts Login Screen will be displayed, then “**Log In**”.

3) You will also see the Patient's information displayed in “**EasyDent Pop Up Window**”. This window will stay on top of the screen for ease of use.

The screenshot shows the Allscripts ePrescribe interface. At the top, the patient's name is "Flintstone, Fred" and the patient ID is "00002". A red arrow labeled "4" points to the "Patient Info 00002" pop-up window, which displays the following information:

| | | | |
|---------------------------|-------------------|------------|------------|
| 00002 | FLINTSTONE FRED N | | |
| Acc Num | Patient Name | | |
| FLINTSTONE | FRED | N | 11/18/1970 |
| Last Name | First Name | In | DOB |
| 23 Walnut Lane | | | |
| SPRINGFIELD | | IL | 62001 |
| City | State | Zip Code | |
| 314 1112222 | 573 2223333 | | |
| Home Phone | Work Phone | Cell Phone | |
| AETNA | | | |
| Primary Insurance Carrier | | | |
| eMail Address | | | |

At the bottom of the main interface, there is a table with the following data:

| Patient ID | Patient Name | DOB | Phone Number | Street Address |
|------------|------------------|------------|---------------|---------------------|
| 02134 | Rubio, Marco | 03/18/1974 | (314)888-5555 | 2900 White House Rd |
| 00002 | Flintstone, Fred | 11/18/1970 | (314)111-2222 | 23 Walnut Lane |

A red arrow labeled "5" points to the row for Fred Flintstone in this table.

4) You can simply click on any field in the “**EasyDent Pop Up Window**”, and the information will automatically be placed in the Windows Clipboard.

5) Next click on any field on the Allscripts screen, right-click and paste.

This lets you easily add your EasyDent patients to the Allscripts system **without re-typing** any patient demographic information.

Once a patient is in Allscripts, you can easily find them by their patient account number or name (which will automatically be displayed in the “**EasyDent Pop Up Window**”).